



Nordson MEDICAL  
805 West 71st Street  
Loveland, CO 80538, USA

+1.888.404.5837  
+1.970.267.5200  
nordsonmedical.com

## **Product Change Notification**

**Product Change #: 1072988**

**Posting Date: 1/11/2019**

**Description of Change: The previously communicated maximum back pressure of 8bar has been redefined at 5 bar, for the following part numbers: SCV67220, SCV67200/230, SCV67220/ML, SCV67230, SCV67230/220, SCV67230/ML, SCV67FL/220, SCV67FL/230, SCV06265, SCV06018. There will be no physical changes made to these parts.**

**Reason: Review of inspection results over the production lifespan has shown that the current acceptable back pressure is close to the maximum load, and 8bar cannot be guaranteed. The maximum back pressure of these part numbers has been redefined at 5bar. However, compliance with ISO 8536-12 (chapter 6.5, backflow of 2bar for 15 minutes) will be maintained**

**Parts Affected: SCV67220, SCV67200/230, SCV67220/ML, SCV67230, SCV67230/220, SCV67230/ML, SCV67FL/220, SCV67FL/230, SCV06265, SCV06018**

**Effective Date: 1/11/2019**

**Implementation Procedure: This change aligns the communicated specifications with actual performance.**