



Nordson MEDICAL
805 West 71st Street
Loveland, CO 80538, USA

+1.888.404.5837
+1.970.267.5200
nordsonmedical.com

Product Change Notification

Date of Correspondence	Planned Date of Change Implementation
<i>3/10/2020</i>	<i>Immediate</i>
Change Notification Number	Sub-Supplier Change Notification Number
<i>1100641 REV A</i>	<i>N/A</i>

Product(s) Affected

Part Number(s)			
<i>FCFM-2A-004</i>			
Change Description			
<i>Product is no longer available.</i>			
Reason for Change			
<i>Through-hole is susceptible to obstruction as the Nylon body material can shave off during assembly of the diverter. The high probability of flaking flash does not allow this product to consistently meet Nordson MEDICAL quality standards.</i>			
Change Level			
Level 1 <input type="checkbox"/>	Level 2 <input type="checkbox"/>	Level 3 <input type="checkbox"/>	Emergency <input checked="" type="checkbox"/>

Implementation Procedure
<i>Product is subject to immediate discontinuation.</i>

Single Point of Contact	Email Address
Customer Service	sales@nordsonmedical.com