



Nordson MEDICAL
805 West 71st Street
Loveland, CO 80538, USA

+1.888.404.5837
+1.970.267.5200
nordsonmedical.com

November 5th, 2021

Product Change Notification - 1307016

Affected Part(s): All products with the suffix -J1A or -J1A2 (Kynar)

Dear Valued Customer:

You are receiving this notice because your company has purchased one or more of the Part Numbers listed above within the last 24 months.

The manufacturer of the Kynar® 1000HD (used in Nordson MEDICAL “-J1A” (natural) or “-J1A2” (black) products) will stop producing this grade by December 31, 2021, and as such, the material will not be available beyond that date. An alternate grade (Kynar® 740) is available as a replacement to the Kynar® 1000HD. The Kynar® 740 polymer is designed to give the same performance characteristics of the Kynar® 1000HD, however, has a slightly darker appearance after processing.

Nordson MEDICAL always strives to give as much notice as possible for changes that impact raw material, however the amount of material sourced is only enough to complete existing orders. **All future orders will need to be for the Kynar® 740 alternate material, indicated with a “-J2A” (natural) or “-J2A2” (black) material designation.**

All processes performed by Nordson MEDICAL will remain unchanged, including the manufacturing environments, production processes, traceability, quality processes, packaging, product specifications and shipment documentation will remain unchanged.

For answers to any questions regarding this change, or if more details are needed, please contact your Strategic Account Manager, our customer service team at (888) 404-5837, or sales@nordsonmedical.com. We appreciate your business and want to thank you for choosing Nordson MEDICAL.

Sincerely,

Kyle R. Steele, P.E.
Global Product Line Manager
Engineered Solutions & Flow Control
Nordson MEDICAL

Affected Products

1050-J1A	FTL220-J1A	M6250-J1A	N680-J1A	T250-J1A
1600-J1A	FTL230-J1A	MTLL210-J1A	N680/655-J1A	T30-J1A
16220-J1A	FTL240-J1A	MTLL220-J1A	PIP070-J1A	T360/230-J1A
16230-J1A	FTLC-J1A	MTLL230-J1A	PIP220-J1A	T40-J1A
16250-J1A	FTLL035-J1A	MTLL240-J1A	PIP230-J1A	T410-J1A
1800-J1A	FTLL210-J1A	MTLL250-J1A	PIP30-J1A	T420-J1A
18210-J1A	FTLL220-J1A	MTLL360-J1A	PIP50-J1A	T420/430-J1A
18220-J1A	FTLL230-J1A	MTLL420-J1A	PIP60-J1A	T430-J1A
18230-J1A	FTLL240-J1A	MTLL430-J1A	PMS230-210-J1A	T430/420-J1A
18420-J1A	FTLL250-J1A	MTLL450-J1A	RS-J1A	T50-J1A
18430-J1A	FTLL360-J1A	MTLLP-J1A	S210-J1A	T60-J1A
1850-J1A	FTLL420-J1A	MTLS210-J1A	S230-J1A	T680-J1A
1860-J1A	FTLL430-J1A	MTLS230-J1A	S250-J1A	X10-J1A
18FTLL-J1A	FTLL445-J1A	MTLS240-J1A	S420-J1A	X210-J1A
18L230-J1A	FTLL450-J1A	MTLSP-J1A	SFM100/SFM075-J1A	X220-J1A
18L260-J1A	FTLLB210-J1A	N004-J1A	SFM150/SFM075-J1A	X230-J1A
18LFTLL-J1A	FTLLB220-J1A	N007/004-J1A	SFM150/SFM100-J1A	X240-J1A
3040-J1A	FTLLB230-J1A	N013-J1A	SFM200/SFM075-J1A	X445-J1A
3050-J1A	FTLLC-J1A	N013/007-J1A	SFM200/SFM150-J1A	XFTL-J1A
3060-J1A	FTLLE-J1A	N025/004-J1A	SFMN6100-J1A	XFTLL-J1A
4040-J1A	FTLLP-J1A	N025/013-J1A	SFMN655-J1A	XL210-J1A
4050-J1A	FTLP-J1A	N035-J1A	SFMN670-J1A	XL220-J1A
4060-J1A	K210-J1A	N035/013-J1A	SFMN680-J1A	XL230-J1A
5060-J1A	K220-J1A	N035/025-J1A	SFMNC-J1A	XL420-J1A
6060-J1A	K230-J1A	N055-J1A	SFMX6100-J1A	XT210-J1A
AA-J1A	K240-J1A	N055/013-J1A	SFMX6110-J1A	Y090-J1A
ABR004-J1A-J1A	K430-J1A	N055/035-J1A	SFMX655-J1A	Y100-J1A
ABR013-J1A-J1A	KL220-J1A	N065-J1A	SFMX670-J1A	Y210-J1A
AC-J1A	KP-J1A	N070/013-J1A	SFMX680-J1A	Y220-J1A
AD-J1A	L10-J1A	N070/055-J1A	SFMXC-J1A	Y230-J1A
BSFTLL-J1A	L210-J1A	N210-J1A	SFTLL-J1A	Y240-J1A
BSMLRL-J1A	L230-J1A	N220/210-J1A	SMTLL-J1A	Y250-J1A
CC-J1A	L230/210-J1A	N230-J1A	SP-J1A	Y260-J1A
CF906-J1A	L230/220-J1A	N230/220-J1A	T10-J1A	Y410-J1A
CF920-J1A	L40-J1A	N410-J1A	T210-J1A	Y420-J1A
CF934-J1A	L420-J1A	N420-J1A	T220-J1A	Y430-J1A
DD-J1A	L440-J1A	N420/410-J1A	T230-J1A	Y450-J1A
FTL10-J1A	LC78-J1A	N430/410-J1A	T230/210-J1A	Y460-J1A
FTL20-J1A	LNS-J1A	N440-J1A	T230/220-J1A	Y630/655-J1A
FTL210-J1A	LP4-J1A	N670-J1A	T230R220-J1A	Y670-J1A
				Y680-J1A



Nordson MEDICAL
805 West 71st Street
Loveland, CO 80538, USA

+1.888.404.5837
+1.970.267.5200
nordsonmedical.com