



**Nordson MEDICAL**  
805 West 71st Street  
Loveland, CO 80538, USA

+1.888.404.5837  
+1.970.267.5200  
nordsonmedical.com

## Product Change Notification

Product Change #: 334

Posting Date: June 29, 2016

1. Description of Change: Change in drawing to match product dimensions. Change the Thru Hole from .117" +/- .005" to .112" +/- .005"

Reason: Have been running under deviation.

Parts Affected: BPC240

Effective Date: Immediately

Implementation Procedure: None. Document change only.