

## Product Change Notification

Product Change #: PCN346

Posting Date: 9-7-2017

Description of Change: A change to the mold relating to the injection subgating to make this part more easily manufactured. This manufacturability change will result in the flow path diameter to change from dia. 0.068" to dia. 0.058" and a cosmetic change to the outward appearance of the part.

Reason: Changes needed to meet production growth.

Parts Affected: N420-1, N420-2, N420-20, N420-J1A2, N420-6005, N420-J1A, N430/420-1, N430/420-2, N430/420-6005, N430/420-6012, N430/420-J1A, N430/410-1, N430/410-2, N430/410-6005, and N430/410-J1A

Effective Date: Production expected by mid-November, with samples projected to be available between late September and early October.

Implementation Procedure: PCN will be sent out in August with a 60-90-day notice period. New products will be produced and move into inventory when the mold improvement is completed. The current product version will ship until inventory is exhausted. New products will then begin shipping, trackable by Lot Number. Note, we will endeavor to ship orders complete with a single lot number.

Part prior to grip change



Part with new subgate detail

