



Nordson MEDICAL
805 West 71st Street
Loveland, CO 80538, USA

+1.888.404.5837
+1.970.267.5200
nordsonmedical.com

Product Change Notification

Product Change #: PCN366

Posting Date: **August 30, 2017**

Description of Change: The ID snap will be called out differently to more closely represent how part is being inspected. Physical part will not change, only how the dimension is called out on the print.


Reason: Parts have been sent to MRB for the two leg's ID being above USL.

Parts Affected: Y6100-6005, Y6110-6005


Effective Date: **September 30, 2017**

Implementation Procedure: Rolling change, no inventory will be scrapped.


Y6100 ID Callout Pre-Change



Y6100 ID Callout Post-Change




Y6110 ID Callout Pre-Change



SECTION A-A

Y6110 ID Callout Post-Change



SECTION A-A

SECTION B-B