


Product Change Notification

Product Change #: PCN382

Posting Date: 05/04/2018

Description of Change:

N430/420 (reducers) parts will have a change to the through hole diameter to enhance flow; this correspondingly results in a slight change to the 420 barb geometry as depicted in the image below.



Reason: Flow enhancement

Parts Affected: N430/420-1, N430/420-2, N430/420-6005,
N430/420-6012, N430/420-J1A

Effective Date: May 2018



Nordson MEDICAL
805 West 71st Street
Loveland, CO 80538, USA

+1.888.404.5837
+1.970.267.5200
nordsonmedical.com

Implementation Procedure: Samples will be available starting May and this product will experience a rolling change with traceability through lot numbers. Old part will be obsoleted after depletion of inventory. Updated engineering drawings will be available on our website.