

**Value Plastics dba Nordson MEDICAL**

805 W 71st Street - Loveland, CO 80538
 Telephone: (970) 267-5200 - Fax: (970)223-0953
 Toll Free in the USA: (888)404-5837
<http://www.nordsonmedical.com>
 Federal Tax ID #:84-0700506

CERTIFICATE OF CONFORMANCE

Material Traceability & Compliance
 Generated on: 5/5/2017

Customer:	EXAMPLE CUSTOMER NAME
Customer NO:	EXAMPLE CUSTOMER NUMBER
PO #:	EXAMPLE PO
Order #:	EXAMPLE ORDER #

Notes:

Value Plastics, INC. dba Nordson MEDICAL, certifies that the products referenced below are manufactured to comply with all technical specifications stated on the Value Plastics website, stated in the current catalog and referenced on current part drawings found at www.valueplastics.com.
 Value Plastics Certifies:

PLEASE NOTE – ALL PARTS WITH LOTS GREATER THAN 150000 WERE MANUFACTURED IN OUR NEW LOVELAND FACILITY. PLEASE SEE THE MASTER VALIDATION PLAN LOCATED ON OUR WEBPAGE http://www.nordsonmedical.com/news/mfg_facility_location.aspx FOR DETAILS RELATED TO THIS CHANGE.

REQUIREMENTS	COMPLIANT
Product Inspections	All products manufactured by Value Plastics are visually inspected to a 0.65 AQL. Value Plastics sample sizes are based on ANSI/ASQCZ 1.4 Table II-A single-sampling plans. Value Plastics does not perform C=0 inspections unless specified.
Manufacturing Process	Value Plastics injection molding processes and secondary operations are conducted in an environment that meets ISO standard 14644-1 Class 8 for clean rooms.
Regulatory Requirements	A comprehensive list of the most up to date regulatory compliance information is available at www.valueplastics.com/tech/material_information.aspx Because of possible changes in the law and/or in regulations, we encourage our customers to periodically verify the status of the regulatory compliance. Please be informed that Value Plastics relies on information provided by our suppliers and vendors, whose materials make up the sole content of our products. There are no processing agents or mold releases used in the manufacture of our products. Value Plastics does not test or analyze these materials for ANY specified regulatory requirements, therefore, customers and end-users must make their own determinations ensuring the use of these products are safe, lawful and suitable for their intended applications.
Change Notification	Most Value Plastics, Inc. products are our proprietary designs. Visit www.valueplastics.com/tech/product_changes_notification.aspx for details of our product notification policy.

Jordan Berkel

Quality Assurance Manager, Nordson MEDICAL

Part #	Description	Customer Part #	Order Lot #	QTY
MTLSP-6005 Other Requirements:	Male Luer Slip Plug, Closed at Grip, Animal-Free Natural Polypropylene	EXAMPLE	EXAMPLE	1000
	Material Lot: EXAMPLE Percent: 100 Material: Flint Hills Resources P5M6K-080 Nat PP - Resin manufacturer has certified the plastic resin to USP Class VI		(6005)	
Part #	Description	Customer Part #	Order Lot #	QTY
MTLL230-J1A Other Requirements:	Male Luer Integral Lock Ring to 200 Series Barb, 1/8" (3.2 mm) ID Tubing, Natural Kynar PVDF	EXAMPLE	EXAMPLE	500
	Material Lot: EXAMPLE Percent: 100 Material: Arkema Kynar (R) 1000 HD Natural PVDF Resin, USP Class VI		(J1A)	



Value Plastics dba Nordson MEDICAL

805 W 71st Street - Loveland, CO 80538
Telephone: (970) 267-5200 - Fax: (970)223-0953
Toll Free in the USA: (888)404-5837
<http://www.nordsonmedical.com>
Federal Tax ID #:84-0700506

CERTIFICATE OF ORIGIN	
Material Traceability & Compliance Generated on: 5/5/2017	
Customer:	EXAMPLE CUSTOMER
Customer NO:	EXAMPLE CUSTOMER NUMBER
PO #:	EXAMPLE PO
Order #:	EXAMPLE ORDER #

The below details and statements are correct, all the goods were produced in their stated country of origin, and they comply with the export administration requirements.

Customer part numbers are given without the guarantee of accuracy.

Tim Faatz

Materials Manager

Part Number	Customer Part Number	Country Of Origin
MTLSP-6005	EXAMPLE	UNITED STATES
MTLL230-J1A	EXAMPLE	UNITED STATES