

**Value Plastics dba Nordson MEDICAL**

805 W 71st Street - Loveland, CO 80538  
 Telephone: (970) 267-5200 - Fax: (970)223-0953  
 Toll Free in the USA: (888)404-5837  
<http://www.nordsonmedical.com>  
 Federal Tax ID #:84-0700506

**CERTIFICATE OF CONFORMANCE**

Material Traceability & Compliance  
 Generated on: 5/5/2017

<b>Customer:</b>	EXAMPLE CUSTOMER NAME
<b>Customer NO:</b>	EXAMPLE CUSTOMER NUMBER
<b>PO #:</b>	EXAMPLE PO
<b>Order #:</b>	EXAMPLE ORDER #

**Notes:**

Value Plastics, INC. dba Nordson MEDICAL, certifies that the products referenced below are manufactured to comply with all technical specifications stated on the Value Plastics website, stated in the current catalog and referenced on current part drawings found at [www.valueplastics.com](http://www.valueplastics.com).  
 Value Plastics Certifies:

PLEASE NOTE – ALL PARTS WITH LOTS GREATER THAN 150000 WERE MANUFACTURED IN OUR NEW LOVELAND FACILITY. PLEASE SEE THE MASTER VALIDATION PLAN LOCATED ON OUR WEBPAGE [http://www.nordsonmedical.com/news/mfg\\_facility\\_location.aspx](http://www.nordsonmedical.com/news/mfg_facility_location.aspx) FOR DETAILS RELATED TO THIS CHANGE.

REQUIREMENTS	COMPLIANT
<b>Product Inspections</b>	All products manufactured by Value Plastics are visually inspected to a 0.65 AQL. Value Plastics sample sizes are based on ANSI/ASQCZ 1.4 Table II-A single-sampling plans. Value Plastics does not perform C=0 inspections unless specified.
<b>Manufacturing Process</b>	Value Plastics injection molding processes and secondary operations are conducted in an environment that meets ISO standard 14644-1 Class 8 for clean rooms.
<b>Regulatory Requirements</b>	A comprehensive list of the most up to date regulatory compliance information is available at <a href="http://www.valueplastics.com/tech/material_information.aspx">www.valueplastics.com/tech/material_information.aspx</a> Because of possible changes in the law and/or in regulations, we encourage our customers to periodically verify the status of the regulatory compliance.  Please be informed that Value Plastics relies on information provided by our suppliers and vendors, whose materials make up the sole content of our products. There are no processing agents or mold releases used in the manufacture of our products. Value Plastics does not test or analyze these materials for ANY specified regulatory requirements, therefore, customers and end-users must make their own determinations ensuring the use of these products are safe, lawful and suitable for their intended applications.
<b>Change Notification</b>	Most Value Plastics, Inc. products are our proprietary designs. Visit <a href="http://www.valueplastics.com/tech/product_changes_notification.aspx">www.valueplastics.com/tech/product_changes_notification.aspx</a> for details of our product notification policy.

**Jordan Berkel**

Quality Assurance Manager, Nordson MEDICAL

Part #	Description	Customer Part #	Order Lot #	QTY
<b>MTLSP-6005</b> Other Requirements:	<b>Male Luer Slip Plug, Closed at Grip, Animal-Free Natural Polypropylene</b>	<b>EXAMPLE</b>	<b>EXAMPLE</b>	<b>1000</b>
	Material Lot: EXAMPLE      Percent: 100	Material: Flint Hills Resources P5M6K-080 Nat PP - Resin manufacturer has certified the plastic resin to USP Class VI	(6005)	
Part #	Description	Customer Part #	Order Lot #	QTY
<b>MTLL230-J1A</b> Other Requirements:	<b>Male Luer Integral Lock Ring to 200 Series Barb, 1/8" (3.2 mm) ID Tubing, Natural Kynar PVDF</b>	<b>EXAMPLE</b>	<b>EXAMPLE</b>	<b>500</b>
	Material Lot: EXAMPLE      Percent: 100	Material: Arkema Kynar (R) 1000 HD Natural PVDF Resin, USP Class VI	(J1A)	



**Value Plastics dba Nordson MEDICAL**

805 W 71st Street - Loveland, CO 80538  
Telephone: (970) 267-5200 - Fax: (970)223-0953  
Toll Free in the USA: (888)404-5837  
http://www.nordsonmedical.com  
Federal Tax ID #:84-0700506

<b>CERTIFICATE OF ORIGIN</b>	
Material Traceability & Compliance Generated on: 5/5/2017	
<b>Customer:</b>	EXAMPLE CUSTOMER
<b>Customer NO:</b>	EXAMPLE CUSTOMER NUMBER
<b>PO #:</b>	EXAMPLE PO
<b>Order #:</b>	EXAMPLE ORDER #

The below details and statements are correct, all the goods were produced in their stated country of origin, and they comply with the export administration requirements.

Customer part numbers are given without the guarantee of accuracy.

**Tim Faatz**

Materials Manager

Part Number	Customer Part Number	Country Of Origin
MTLSP-6005	EXAMPLE	UNITED STATES
MTLL230-J1A	EXAMPLE	UNITED STATES