



MQC SERIES™

QUICK CONNECT FITTINGS

Nordson MEDICAL's MQC Series™ of quick connects not only have an intuitive push to connect design, but the variety of colors and barb sizes allows you to create your own custom solution. Between the large selection of colors and the MQC's unique, innovative design, misconnections will be a problem of the past!



FEATURES

Push-to-Connect design

Large, ergonomic buttons

Audible click on connect

Lubricated O-ring seal

Unique design & variety of colors

BENEFITS

Intuitive and quick connection

Easy to use with gloved hands

Confidence in a secure connection

Long service life

Helps to prevent misconnections

SPECIFICATIONS

Materials: Standard Nylon, Buna-N; USP Class VI Nylon, USP Class VI Silicone

Sterilization: Gamma radiation stable, EtO compatible

Barb Size: 200 Series - 1/8" (3.2 mm), 5/32" (4.0 mm), 3/16" (4.8 mm)

Colors: White, Red, Blue, Green, Purple, Black, Grey

Leak Resistance: < .007 psi/0.05 kPa/sec

Tensile Pull Load: 25 lbs

Life Cycles: Up to 1,500


Working Pressure: Up to 30 psi

For more information and samples, please visit www.nordsonmedical.com


MQCs have ergonomic buttons for an easy connect and disconnect every time. Listen for the loud "click" and a secure connection will never be doubted. Contact us today for a free sample!

MQC SERIES™


QUICK CONNECT FITTINGS




MQCRM230-001
Dual Latch Quick Connect Male
200 Series Barb
1/8" (3.2 mm) ID Tubing, White Nylon
Dimensions A: 1.000" B: 2.011" C: .528"
Also available in colors listed below




MQCF230-1
Dual Latch Quick Connect Female
200 Series Barb
1/8" (3.2 mm) ID Tubing, White Nylon
A: 1.000" B: 2.011" C: .493"
Also available in colors listed below




MQCRM240-001
Dual Latch Quick Connect Male
200 Series Barb
5/32" (4.0 mm) ID Tubing, White Nylon
A: 1.000" B: 2.203" C: .528"
Also available in colors listed below



MQCF240-1
Dual Latch Quick Connect Female
200 Series Barb
5/32" (4.0 mm) ID Tubing, White Nylon
A: 1.000" B: 2.203" C: .493"
Also available in colors listed below





MQCRM250-001
Dual Latch Quick Connect Male
200 Series Barb
3/16" (4.8 mm) ID Tubing, White Nylon
A: 1.000" B: 2.399" C: .528"
Also available in colors listed below



MQCF250-1
Dual Latch Quick Connect Female
200 Series Barb
3/16" (4.8 mm) ID Tubing, White Nylon
A: 1.000" B: 2.399" C: .493"
Also available in colors listed below

Flow Rate vs. Pressure (Water)
MQC Series™ 240 Barb
Upstream Pressure (Bar)



Color Availability - Standard Nylon



Color Availability - USP Class VI Nylon



The data presented here are for reference only. They are not presented as a specification guideline. As your application conditions are beyond Nordson MEDICAL's control, it is imperative that you test products in their specific environments to determine suitability for a particular use. For further information on Nordson MEDICAL warranties, please review section 10 of the "Terms and Conditions of Sale" on the Nordson MEDICAL website at http://www.nordsonmedical.com/technical/orders/vp_terms.pdf or Nordson MEDICAL's most recent catalog. Patent pending.

United States of America

3325 South Timberline Road
Fort Collins, CO 80525
U.S.A.

Phone: +1 (970) 267.5200
Toll Free: (888) 404.5837
sales@nordsonmedical.com
www.nordsonmedical.com

Asia Pacific

Rm. 206, 2/F, Tung Ming Building,
40-42 Des Voeux Road Central
Hong Kong

Phone: (852) 6808 9862
Fax: (852) 2470 9802
cnsales@nordsonmedical.com
www.nordsonmedical.com

IN-STOCK

Components & Tools



Experts in Medical Device Design, Manufacturing & Component Technologies


Browse our online store, with more than 7,000 stock components and technologies, most ready to ship within 24 hours.




Components & Technologies




Catheters & Cannulae




Medical Balloons



Medical Tubing



Fluid Management Components



Biomaterial Delivery Devices



Nitinol Components