

Compliant to the ISO 80369-6 Standard

NE Series connectors are designed, tested, qualified, and documented to be compliant with the soon-to-be published ISO 80369-6 standard. Small and lightweight, these connectors are easy-to-clean and offer comfort and compliance. The medical-grade material is clear with no colorants added to simplify regulatory management and also features excellent solvent bonding properties. The female connectors showcase wings and the male connectors have an easy-to-grip hub. The available caps and plugs seal off the connectors for irradiation, and there are vented versions for use with ethylene oxide (EtO) sterilization.



Features & Benefits


- ISO 80369-6 compliant
- Small, lightweight design
- Caps available vented and non-vented
- Easy to use bond-in tube pockets
- Regulatory information available
- Clear materials with no colorants
- Custom packaging available
- Ergonomic wings on female connectors
- Male connectors with an easy-to-grip hub

Specifications


- **Materials:**
Medical-grade acrylic resin
Custom colorants and materials available
- **Sizes:**
Male and female with 2.1 mm or 3.1 mm pocket sizes
Custom pocket sizes available
- **Assembly Sets:**
Readily available for custom options
- **Sterilization:**
Gamma radiation stable, EtO compatible

NE Series


Neuraxial Connectors




NEFE21 (-50)
Neuraxial Female Connector with 2.1 mm Tubing
Bond-in Pocket in Natural Acrylic




NEME21 (-50)
Neuraxial Male Connector with 2.1 mm Tubing
Bond-in Pocket in Natural Acrylic




NEFE31 (-50)
Neuraxial Female Connector with 3.1 mm Tubing
Bond-in Pocket in Natural Acrylic



NEME31 (-50)
Neuraxial Male Connector with 3.1 mm Tubing
Bond-in Pocket in Natural Acrylic




NEFC (-50)
Neuraxial Female Cap in Natural Acrylic




NEMP (-50)
Neuraxial Male Plug in Natural Acrylic

ISO 80369-6 Accessories

Vented caps and plugs are not covered under the ISO 80369-6 standard, but are compatible with our neuraxial male and female connectors.



NEFCV (-50)
Neuraxial Female Vented Cap in Natural Acrylic



NEMPV (-50)
Neuraxial Male Vented Plug in Natural Acrylic

MAR-NE-DS-01

IN-STOCK

Components & Tools



Experts in Medical Device Design, Manufacturing & Component Technologies


Browse our online store, with more than 7,000 stock components and technologies, most ready to ship within 24 hours.



Components & Technologies




Catheters & Cannulae




Medical Balloons



Medical Tubing



Fluid Management Components



Biomaterial Delivery Devices



Nitinol Components