



The XQ Series quick connect fittings/quick disconnect couplings, bring performance, convenience and utility to your tubing applications. Offered in three tubing sizes to fit 1/4" (6.4 mm), 5/16" (8.0 mm) and 3/8" (9.5 mm) ID flexible tubing, XQ series products may be used with PVC, polyurethane, silicone and other types of flexible tubing. Nordson MEDICAL's exceptional single barb designs provide superior pulloff resistance when compared to multi-barb styles. For versatility, durable XQ Series couplings interchange with other plastic and metal button-type couplings. XQ Series quick connects are a perfect fit for your tubing application.

Features

Integral automatic shutoff valves
Buttons and pin springs
Audible connection signal
Ergonomic grip
700 Series single barb design
Fits 1/4" (6.4 mm), 5/16" (8.0 mm)
and 3/8" (9.5 mm) ID flexible tubing

Materials

Acetal copolymer resin
Polypropylene resin
Buna-N o-ring elastomer seals
EPDM o-ring elastomer seals
Stainless steel valve springs

Specifications

Maximum pressure:
120 psi, 8.28 bar

Leak resistance:
< .007 psi / 0.05 kPa / sec.

Maximum C_v:
0.58 - 1/4" (6.4 mm), both sides shutoff
.76 - 3/8" (9.5 mm), both sides shutoff

Maximum Temperature:
180°F, (82.2°C) to -40°F, (-40°C)

For more information and samples, please visit www.nordsonmedical.com

Part Number

Description

All valved connectors are 100% leak tested!

XQSVF755-1006-B	Female Valved Connector, 1/4" (6.4 mm) ID Tubing, Natural Acetal	
XQSVF765-1006-B	Female Valved Connector, 5/16" (8.0 mm) ID Tubing, Natural Acetal	
XQSVF770-1006-B	Female Valved Connector, 3/8" (9.5 mm) ID Tubing, Natural Acetal	
XQSVF755-1007-B	Female Valved Connector, 1/4" (6.4 mm) ID Tubing, Black Acetal	
XQSVF765-1007-B	Female Valved Connector, 5/16" (8.0 mm) ID Tubing, Black Acetal	
XQSVF770-1007-B	Female Valved Connector, 3/8" (9.5 mm) ID Tubing, Black Acetal	
XQSVF755-6006-EP	Female Valved Connector, 1/4" (6.4 mm) ID Tubing, White Polypropylene	
XQSVF770-6006-EP	Female Valved Connector, 3/8" (9.5 mm) ID Tubing, White Polypropylene	
XQSF755-1006	Female Open Flow Connector, 1/4" (6.4 mm) ID Tubing, Natural Acetal	
XQSF770-1006	Female Open Flow Connector, 3/8" (9.5mm) ID Tubing, Natural Acetal	
XQSF770-1007	Female Open Flow Connector, 700 Series Barbs, 3/8" (9.5mm) ID Tubing, Black Acetal	
XQSVBF755-1006-B	Female Panel Mount Valved Connector, 1/4" (6.4 mm) ID Tubing, Natural Acetal	
XQSVBF770-1006-B	Female Panel Mount Valved Connector, 3/8" (9.5 mm) ID Tubing, Natural Acetal	
XQSBF755-1006	Female Panel Mount Open Flow Connector, 1/4" (6.4 mm) ID Tubing, Natural Acetal	
XQSBF770-1006	Female Panel Mount Open Flow Connector, 3/8" (9.5 mm) ID Tubing, Natural Acetal	
XQCVM755-1006-B	Male Valved Connector, 1/4" (6.4 mm) ID Tubing, Natural Acetal	
XQCVM765-1006-B	Male Valved Connector, 700 Series Barb, 5/16" (8.0 mm) ID Tubing, Natural Acetal	
XQCVM770-1006-B	Male Valved Connector, 3/8" (9.5 mm) ID Tubing, Natural Acetal	
XQCVM765-1007-B	Male Valved Connector, 700 Series Barb, 5/16" (8.0 mm) ID Tubing, Black Acetal	
XQCVM770-1007-B	Male Valved Connector, 700 Series Barb, 3/8" (9.5 mm) ID Tubing, Black Acetal	
XQCVM755-6006-EP	Male Valved Connector, 1/4" (6.4 mm) ID Tubing, White Polypropylene	
XQCVM770-6006-EP	Male Valved Connector, 3/8" (9.5 mm) ID Tubing, White Polypropylene	
XQCM755-1006-B	Male Open Flow Connector, 1/4" (6.4 mm) ID Tubing, Natural Acetal	
XQCM770-1006-B	Male Open Flow Connector, 3/8" (9.5 mm) ID Tubing, Natural Acetal	
XQCVBM755-1006-B	Male Panel Mount Valved Connector, 1/4" (6.4 mm) ID Tubing, Natural Acetal	
XQCVBM770-1006-B	Male Panel Mount Valved Connector, 3/8" (9.5 mm) ID Tubing, Natural Acetal	
XQCBM755-1006-B	Male Panel Mount Open Flow Connector, 1/4" (6.4 mm) ID Tubing, Natural Acetal	
XQCBM770-1006-B	Male Panel Mount Open Flow Connector, 3/8" (9.5 mm) ID Tubing, Natural Acetal	
XQCVLM755-1006-B	Elbow Male Valved Connector, 1/4" (6.4 mm) ID Tubing, Natural Acetal	
XQCLM755-1006-B	Elbow Male Open Flow Connector, 1/4" (6.4 mm) ID Tubing, Natural Acetal	
XQCLM770-1006-B	Elbow Male Open Flow Connector, 3/8" (9.5 mm) ID Tubing, Natural Acetal	

B = Buna-N o-ring

EP = EPDM o-ring

IN-STOCK

Components & Tools



Experts in Medical Device Design, Manufacturing & Component Technologies


Browse our online store, with more than 7,000 stock components and technologies, most ready to ship within 24 hours.




Components & Technologies




Catheters & Cannulae




Medical Balloons



Medical Tubing



Fluid Management Components



Biomaterial Delivery Devices



Nitinol Components