




NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE
DIMENSIONS MARKED \varnothing ARE CONTROL DIMENSIONS

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	INITIAL RELEASE	06-12-01	
B	CO 015 Drafting Update	07-11-05	Rick C.
C	ECO 227: Obsolete Profax PD626 (-6) Polypropylene	11/9/12	S.N.R
D	ECO 285	11-13-14	H.L

ITEM NO.	PART NUMBER	DESCRIPTION	QTY.
1	FLLR	Finger Snap Luer Lock Ring	1
2	BSMLRL	Male Luer with 5/16" Hex to 1/4-28 UNF Bottom Sealing Thread	1



SECTION A-A

Material Identification			
Assembly Number	Item 1 (FLLR) Material	Item 2 (BSMLRL) Material	Description
ABSML-20-06	Orange Nylon	Kynar PVDF	ABSML-Orange Nylon Over Kynar PVDF
ABSML-20-07	Yellow Nylon	Kynar PVDF	ABSML-Yellow Nylon Over Kynar PVDF
ABSML-20-1	White Nylon	Kynar PVDF	ABSML-White Nylon Over Kynar PVDF
ABSML-20-2	Black Nylon	Kynar PVDF	ABSML-Black Nylon Over Kynar PVDF
ABSML-20-3	Red Nylon	Kynar PVDF	ABSML-Red Nylon Over Kynar PVDF
ABSML-20-4	Green Nylon	Kynar PVDF	ABSML-Green Nylon Over Kynar PVDF
ABSML-20-5	Blue Nylon	Kynar PVDF	ABSML-Blue Nylon Over Kynar PVDF
ABSML-6-06	Orange Nylon	Polypropylene	ABSML-Orange Nylon Over Polypropylene
ABSML-6-07	Yellow Nylon	Polypropylene	ABSML-Yellow Nylon Over Polypropylene
ABSML-6-1	White Nylon	Polypropylene	ABSML-White Nylon Over Polypropylene
ABSML-6-2	Black Nylon	Polypropylene	ABSML-Black Nylon Over Polypropylene
ABSML-6-3	Red Nylon	Polypropylene	ABSML-Red Nylon Over Polypropylene
ABSML-6-4	Green Nylon	Polypropylene	ABSML-Green Nylon Over Polypropylene
ABSML-6-5	Blue Nylon	Polypropylene	ABSML-Blue Nylon Over Polypropylene
ABSML-6005-06	Orange Nylon	Polypropylene	ABSML-Orange Nylon Over Animal Free PP
ABSML-6005-07	Yellow Nylon	Polypropylene	ABSML-Yellow Nylon Over Animal Free PP
ABSML-6005-1	White Nylon	Polypropylene	ABSML-White Nylon Over Animal Free PP
ABSML-6005-2	Black Nylon	Polypropylene	ABSML-Black Nylon Over Animal Free PP
ABSML-6005-3	Red Nylon	Polypropylene	ABSML-Red Nylon Over Animal Free PP
ABSML-6005-4	Green Nylon	Polypropylene	ABSML-Green Nylon Over Animal Free PP
ABSML-6005-5	Blue Nylon	Polypropylene	ABSML-Blue Nylon Over Animal Free PP
ABSML-6012-2	Black Nylon	Polypropylene	ABSML-Black Nylon Over Black Animal Free PP
ABSML-J1A-06	Orange Nylon	Kynar PVDF	ABSML-Orange Nylon Over Kynar PVDF
ABSML-J1A-07	Yellow Nylon	Kynar PVDF	ABSML-Yellow Nylon Over Kynar PVDF
ABSML-J1A-1	White Nylon	Kynar PVDF	ABSML-White Nylon Over Kynar PVDF
ABSML-J1A-	Black Nylon	Kynar PVDF	ABSML-Black Nylon Over Kynar PVDF
ABSML-J1A-3	Red Nylon	Kynar PVDF	ABSML-Red Nylon Over Kynar PVDF
ABSML-J1A-4	Green Nylon	Kynar PVDF	ABSML-Green Nylon Over Kynar PVDF

Dimension	Material Specific Tolerances							
	Acetal	Nylon	PP	ABS	PC	PVC	Polysulfone	Kynar
A	.312	.312	.312	.312	.312	.312	.312	.307
B								
C								
D								
E								
F								

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES
TOLERANCES ARE:

FRACTIONS	DECIMAL	ANGULAR			
±	.XX ±.01	± 1°			
	.XXX ±.005				
	.XXXX ±.001				

DO NOT SCALE DRAWING

VALUE PLASTICS, INC

PART NUMBER: ABSML

DESCRIPTION: Male Luer with Rotating Lock Ring and 5/16" Hex to 1/4-28 UNF Bottom Sealing Thread

DRAWN BY: Frank.Lombardi	DATE: 11/14/2014	SIZE: B	SCALE: NTS	VER:	DWG. NO. 1006916	REV D
-----------------------------	---------------------	-------------------	---------------	------	---------------------	----------

MANAGER:
Lou Huadong

© 2011 VALUE PLASTICS, INC

SHEET: 1 OF 1

FOR REFERENCE ONLY