



Tubing Specifications

Polyurethane

Part Number	Length	Inside Diameter	Wall Thickness (Ref. Only)	Outside Diameter
UR00-1031C	100', +5', -0'	.062" ± .005"	.031"	.125" ± .005"
UR00-1031M	1000', +5', -0'	.062" ± .005"	.031"	.125" ± .005"
UR00-2047C	100', +5', -0'	.093" ± .005"	.047"	.187" ± .005"
UR00-2047M	1000', +5', -0'	.093" ± .005"	.047"	.187" ± .005"
UR00-3062C	100', +5', -0'	.125" ± .005"	.062"	.250" ± .005"
UR00-3062M	1000', +5', -0'	.125" ± .005"	.062"	.250" ± .005"

LLDPE - Linear Low Density Polyethylene

Part Number	Length	Inside Diameter	Wall Thickness (Ref. Only)	Outside Diameter
LLDPE02-1019C	100', +5', -0'	.062" ± .004"	.019"	.100" ± .004"
LLDPE02-1019M	1000', +5', -0'	.062" ± .004"	.019"	.100" ± .004"
LLDPE02-1031C	100', +5', -0'	.062" ± .004"	.031"	.125" ± .004"
LLDPE02-1031M	1000', +5', -0'	.062" ± .004"	.031"	.125" ± .004"
LLDPE02-2031C	100', +5', -0'	.093" ± .004"	.031"	.156" ± .004"
LLDPE02-2031M	1000', +5', -0'	.093" ± .004"	.031"	.156" ± .004"
LLDPE02-3031C	100', +5', -0'	.125" ± .004"	.031"	.187" ± .004"
LLDPE02-3031M	1000', +5', -0'	.125" ± .004"	.031"	.187" ± .004"

TPE - Thermoplastic Elastomer

Part Number	Length	Inside Diameter	Wall Thickness (Ref. Only)	Outside Diameter
AF7500077	50' ± 1'6"	.125" ± .005"	.062"	.250" ± .010"
AF7500077D	500' ± 15"	.125" ± .005"	.063"	.250" ± .010"
AF7500155	50' ± 1'6"	.188" ± .008"	.063"	.313" ± .011"
AF7500285C	100' ± 1'6"	.250" ± .008"	.063"	.375" ± .011"
AF7500311	50' ± 1'6"	.250" ± .008"	.094"	.438" ± .013"
AF7500506	50' ± 1'6"	.375" ± .011"	.125"	.625" ± .018"
AF7500584	50' ± 1'6"	.500" ± .015"	.125"	.750" ± .020"
AF7500753	50' ± 1'6"	.750" ± .023"	.125"	1.000" ± .030"



Tubing Specifications

Platinum-Cured Silicone

Part Number	Length	Inside Diameter	Wall Thickness (Ref. Only)	Outside Diameter
SILC7001098	50'± 1'6"	.063" ± .005"	.031"	.125" ± .007"
SILC7001280	50'± 1'6"	.078" ± .005"	.031"	.141" ± .007"
SILC7001422	500'± 15"	.063" ± .005"	.063"	.188" ± .010"
SILC7001424	500'± 15"	.125" ± .007"	.063"	.250" ± .010"
SILC7001657	50'± 1'6"	.125" ± .007"	.063"	.250" ± .010"
SILC7001839	100'± 1'6"	.156" ± .007"	.063"	.281" ± .010"
SILC7002203	50'± 1'6"	.250" ± .010"	.063"	.375" ± .015"
SILC7002593	50'± 1'6"	.375" ± .015"	.125"	.625" ± .020"
SILC7002775	50'± 1'6"	.500" ± .020"	.125"	.750" ± .025"
SILC7003139	50'± 1'6"	.750" ± .020"	.125"	1.000" ± .025"

Pump Grade Platinum-Cured Silicone

Part Number	Length	Inside Diameter	Wall Thickness	Outside Diameter (Ref. Only)
SILPG7021044	500'± 15"	.063" ± .005"	.063" ± .004"	.188"
SILPG7021041	300'± 9"	.250" ± .006"	.094" ± .005"	.438"
SILPG7021042	250'± 7.5"	.188" ± .006"	.094" ± .005"	.375"
SILPG7021043	500'± 15"	.125" ± .006"	.063" ± .004"	.250"

PVC - Anti-Microbial

Part Number	Length	Inside Diameter	Wall Thickness (Ref. Only)	Outside Diameter
PVCA1190035	100'± 1'	.188" ± .008"	.063"	.313" ± .010"
PVCA1190049	100'± 1'	.188" ± .008"	.125"	.438" ± .010"
PVCA1190112	100'± 1'	.250" ± .008"	.063"	.375" ± .010"
PVCA1190126	100'± 1'	.250" ± .008"	.125"	.500" ± .010"
PVCA1190189	100'± 1'	.313" ± .008"	.063"	.438" ± .010"
PVCA1190203	100'± 1'	.313" ± .008"	.125"	.563" ± .010"
PVCA1190266	100'± 1'	.375" ± .008"	.063"	.500" ± .010"
PVCA1190280	100'± 1'	.375" ± .008"	.125"	.625" ± .010"
PVCA1190343	100'± 1'	.500" ± .008"	.063"	.625" ± .010"
PVCA1190357	100'± 1'	.500" ± .008"	.125"	.750" ± .010"