



NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP NEW PRODUCT DEVELOPMENT
DIMENSIONS MARKED ARE CONTROL DIMENSIONS
VALUE PLASTICS CPAX300 CLAMPS ARE DESIGNED TO COMPLY WITH ASME BPE-2022
TABLE DT-9.3.1 HYGENIC CLAMP FERRULE: DESIGN CRITERIA

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	INITIAL RELEASE	01-31-23	FJL



PACKAGING POSITION



FOR REFERENCE ONLY

ITEM NO.	PART NUMBER	DESCRIPTION	QTY.
1	CPAX300U	EASY-OPEN CLAMP, 3" (76.2MM), UPPER	1
2	CPAX300L	EASY-OPEN CLAMP, 3" (76.2MM), LOWER	1

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE: FRACTIONS DECIMALS ANGULAR ± .XX ±.01 ±1° .XXX ±.005 .XXXX ±.001 DO NOT SCALE DRAWING			VALUE PLASTICS, INC					
			PART NUMBER: CPAX300					
			DESCRIPTION: TRI-CLOVER EASY-OPEN CLAMP FOR 3" (76.2MM) FLANGE					
DRAWN BY: Jack.Wise		DATE: 2/2/2023		SIZE: C	SCALE: NTS	DRAWING NUMBER: 1312339		REV A
ENG MGR: Frank.Lombardi		DATE: 2/2/2023		© 2014 VALUE PLASTICS, INC				SHEET: SHEET 1 OF 1