



NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP NEW PRODUCT DEVELOPMENT
DIMENSIONS MARKED  ARE CONTROL DIMENSIONS
VALUE PLASTICS CPAX400 CLAMPS ARE DESIGNED TO COMPLY WITH ASME BPE-2019
TABLE DT-9.3-1 HYGENIC CLAMP FERRULE: DESIGN CRITERIA

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	INITIAL RELEASE	02-08-23	



PACKAGING POSITION



FOR REFERENCE ONLY

ITEM NO.	PART NUMBER	DESCRIPTION	QTY.
1	CPAX400H	EASY-OPEN CLAMP, 4" (101.6MM), GENDERLESS HINGE	1
2	CPAX400X	EASY-OPEN CLAMP, 4" (101.6MM), GENDERLESS LATCH	2

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE: FRACTIONS DECIMAL ANGULAR ± ± ± .XX ±.01 .XXX ±.005 .XXXX ±.001 DO NOT SCALE DRAWING		VALUE PLASTICS, INC			
		PART NUMBER: CPAX400			
DESCRIPTION: EASY-OPEN CLAMP, 4" (101.6MM)		SIZE: B	SCALE: NTS	VER: 	DWG. NO. 1299279
DRAWN BY: Jack.Wise		DATE: 2/8/2023		REV A	
MANAGER: Frank.Lombardi		© 2014 VALUE PLASTICS, INC			SHEET: 1 OF 1