





NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE
 DIMENSIONS MARKED \varnothing ARE CONTROL DIMENSIONS

PARTS ARE ASSEMBLED IN THE ORIENTATION SHOWN ON THE DRAWING

ANGLE ALLOWANCE OF DIVERTER ± 10 DEGREES

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	INITIAL RELEASE	08/03/15	K.R.S.
B	ECO 311	01/27/16	K.R.S.
C	ECO 340	07-29-16	FJL



SECTION A-A

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE:		VALUE PLASTICS, INC			
FRACTIONS DECIMAL ANGULAR \pm .XX \pm .01 $\pm 1^\circ$.XXX \pm .005 .XXXX \pm .001		PART NUMBER: FCHPFMR-002			
DO NOT SCALE DRAWING		DESCRIPTION: STOPCOCK, 1-WAY, FL TO ML, LUER LOCK NUT, HIGH PRESSURE			
DRAWN BY: kyles	DATE: 7/29/2016	SIZE: B	SCALE: NTS	VER: C	DWG. NO. 1032423
MANAGER: Frank.Lombardi		© 2014 VALUE PLASTICS, INC			SHEET: 1 OF 1





FOR REFERENCE ONLY

NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE
 DIMENSIONS MARKED $\text{\textcircled{\scriptsize \square}}$ ARE CONTROL DIMENSIONS

$\text{\textcircled{\scriptsize \square}}$ PARTS ARE ASSEMBLED IN THE ORIENTATION SHOWN ON THE DRAWING

$\text{\textcircled{\scriptsize \square}}$ ANGLE ALLOWANCE OF DIVERTER ± 10 DEGREES

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	INITIAL RELEASE	09/16/15	K.R.S.
B	ECO 311	01/27/16	K.R.S.
C	ECO 340	07-29-16	FJL



SECTION A-A

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE:		VALUE PLASTICS, INC			
FRACTIONS DECIMAL ANGULAR \pm .XX \pm .01 $\pm 1^\circ$.XXX \pm .005 .XXXX \pm .001		PART NUMBER: FCHPFMR-003			
DO NOT SCALE DRAWING		DESCRIPTION: STOPCOCK, 1-WAY, FL TO ML, LUER LOCK RING, HIGH PRESSURE			
DRAWN BY: kyles	DATE: 7/29/2016	SIZE: B	SCALE: NTS	VER: 	DWG. NO. REV 1032425 C
MANAGER: Frank.Lombardi		© 2014 VALUE PLASTICS, INC			SHEET: 1 OF 1




FOR REFERENCE ONLY

NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE
 DIMENSIONS MARKED \square ARE CONTROL DIMENSIONS

PARTS ARE ASSEMBLED IN THE ORIENTATION SHOWN ON THE DRAWING

ANGLE ALLOWANCE OF DIVERTER ± 10 DEGREES

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	INITIAL RELEASE	08/19/15	K.R.S.
B	ECO 311	01/27/16	K.R.S.
C	ECO 340	07-29-16	FJL



SECTION A-A

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE:		VALUE PLASTICS, INC			
FRACTIONS	DECIMAL				
±	.XX ±.01	± 1°	PART NUMBER: FCHPFMR-005		
	.XXX ±.005		DESCRIPTION: STOPCOCK, 1-WAY, FL TO ML, LUER LOCK RING, HIGH PRESSURE		
	.XXXX ±.001		DO NOT SCALE DRAWING		
DRAWN BY: kyles	DATE: 7/29/2016	SIZE: B	SCALE: NTS	VER: C	DWG. NO. 1032426
MANAGER: Frank.Lombardi		© 2014 VALUE PLASTICS, INC			SHEET: 1 OF 1

FOR REFERENCE ONLY