



NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE
DIMENSIONS MARKED ARE CONTROL DIMENSIONS

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	INITIAL RELEASE	07/26/04	
B	CO 046 DRAFTING UPDATE	03/17/06	Rick C.
C	ECO 356	06-19-17	FJL



UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE: FRACTIONS DECIMAL ANGULAR ± .XX ±.01 ± 1° .XXX ±.005 .XXXX ±.001 DO NOT SCALE DRAWING		VALUE PLASTICS, INC			
		PART NUMBER: MLBC7			
		DESCRIPTION: LARGE BORE COUPLER W/ INTEGRAL LOCK RING TO MALE SLIP LUER W/ ROTATING LOCK RING (FSLLR)			
DRAWN BY: marie.morleyjahn	DATE: 6/19/2017	SIZE: B	SCALE: NTS	VER: 	DWG. NO. 1005049
NEW PRODUCT DEVELOPMENT MGR: Frank.Lombardi		© 2017 VALUE PLASTICS, INC			REV C
				SHEET: 1 OF 1	

FOR REFERENCE ONLY