




NOTE: FOR MORE DETAILS SEE YOUR V.P.CUSTOMER SERVICE REP.

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	ECO-0029	08-16-00	
B	CO 015 Drafting Update	05-16-05	Rick C.



UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE:		<b>VALUE PLASTICS, INC.</b>	
FRACTIONS ±	DECIMALS .XX ±.01 .XXX ±.005 .XXXX ±.001	ANGULAR ±1°	PART NUMBER: <b>MTLL240</b>
DO NOT SCALE DRAWING		DESCRIPTION: Male Luer Integral Lock Ring to 200 Series Barb, 5/32" (4 mm) ID Tubing	
DRAWN BY: <b>Frank Lombardi 4-18-05</b>	DATE:	SIZE: <b>B</b>	SCALE: <b>NTS</b>
ENGINEER MANAGER: <b>Rick Cairns</b>		DRAWING NUMBER: <b>1005815</b>	REV. <b>B</b>
© 2004 VALUE PLASTICS, INC.			SHEET: <b>1 OF 1</b>

**FOR REFERENCE ONLY**