




NOTE: FOR MORE DETAILS SEE YOUR V.P.CUSTOMER SERVICE REP.

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	ECO-0029	08-18-00	
B	CO 015 Drafting Update	05-16-05	Rick C.



UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE:			VALUE PLASTICS, INC.		
FRACTIONS DECIMALS ANGULAR ± .XX ±.01 ±1° .XXX ±.005 .XXXX ±.001			PART NUMBER: MTLL420		
DO NOT SCALE DRAWING			DESCRIPTION: Male Luer Integral Lock Ring to 400 Series Barb, 3/32" (2.25 mm) ID Tubing		
DRAWN BY: Frank Lombardi 04-27-05			SIZE: B	SCALE: NTS	DRAWING NUMBER: 1006418
Engineer Manager: Rick Cairns					REV. B
© 2005 VALUE PLASTICS, INC.					SHEET: 1 OF 1

FOR REFERENCE ONLY