





NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE
DIMENSIONS MARKED  ARE CONTROL DIMENSIONS

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	INITIAL RELEASE	02-19-97	
B	CO 039 DRAFTING UPDATE	05-08-06	RICK C.
C	ECO 182 DRAFTING UPDATE	04-11-11	RICK C.
D	ECO 346	12-15-17	FJL



UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE: FRACTIONS DECIMAL ANGULAR ± .XX ±.01 ± 1° .XXX ±.005 .XXXX ±.001 DO NOT SCALE DRAWING		VALUE PLASTICS, INC				
		PART NUMBER: N430-410				
		DESCRIPTION: Straight Through Reduction Connector with 400 Series Barbs, 1/8" (3mm) and 1/16" (1.5mm) ID Tubing				
DRAWN BY: kyles	DATE: 08/18/2016	SIZE: B	SCALE: NTS	VER: 	DWG. NO. 1008189	REV D
MANAGER: Frank.Lombardi		© 2011 VALUE PLASTICS, INC				SHEET: 1 OF 1

FOR REFERENCE ONLY