


NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE





SECTION B-B

REVISIONS								
REV.	DESCRIPTION	DATE	APPROVED					
A	INITIAL RELEASE	01/20/09	Rick C.					



UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE:		VALUE PLASTICS, INC							
FRACTIONS DECIMAL ANGULAR  2 .XX 2.01 2 1°  .XXX 2.005  .XXXX 2.001  DO NOT SCALE DRAWING		PART NUMBER:  SCV23050  DESCRIPTION: Check Valve, Female Locking Luer to Male Locking Luer, Clear SAN w/Silicone Diaphragm							
DRAWN BY: MICHELE SCHENK ENGINEERING MANAGER:	DATE: 1/20/2009	B SIZE:	SCALE: NTS	VER:	DWG. NO.	,			REV A
RICK CAIRNS		0	2009	VALUI	E PLASTICS,	INC	SHEET:	1 OF 1	•

## FOR REFERENCE ONLY