
NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE DIMENSIONS MARKED ARE CONTROL DIMENSIONS

REVISIONS								
REV.	DESCRIPTION	DATE	APPROVED					
А	INITIAL RELEASE	06/07/13	S.N.R					





	SCV67220 Bill Of Materials								
POS.	QTY	DESCRIPTION	MATERIAL, COLOR						
1	1	Check Valve V2000 Inlet with Barb 3/32" (2.4 mm) I.D. Tubing	MABS White						
2	1	Check Valve V2000 Outlet with Barb 3/32" (2.4 mm) I.D. Tubing	MABS, Transparent						
3	1	Valve Membrane V2000 IV, 60A Durometer	Silicone, Transparent						
4	рс	Solvent	Tetrahydrofurane						

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE:		VALUE PLASTICS, INC									
FRACTIONS DECIMAL ANGULAR ± .XX ±.01 ± 1° .XXX ±.005 .XXXX ±.005		PART NUMBER: SCV67220									
DO NOT SCALE DRAWING		DESCRIPTION: Check Valve, Barb 3/32" (2.4 mm) ID Tubing to Barb 3/32" (2.4 mm) ID Tubing									
DRAWN BY:	DATE:	SIZE:	SCALE:	VER:	DWG. NO.						REV
kyles	6/7/2013	В	NTS		102717	4					A
MANAGER:											
RaviN		0	2011	VALU	E PLASTICS,	INC	SHEET:	1	OF	1	

FOR REFERENCE ONLY