

NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE DIMENSIONS MARKED ARE CONTROL DIMENSIONS

REVISIONS					
REV.	DESCRIPTION	DATE	APPROVED		
А	INITIAL RELEASE	06/07/13	S.N.R		





SCV67230/ML Bill Of Materials					
POS.	QTY	DESCRIPTION	MATERIAL, COLOR		
1	1	Check Valve V2000 Inlet with Barb 1/8" (3.2 mm) I.D. Tubing	MABS White		
2	1	Check Valve V2000 Outlet with Male Luer Lock	SAN, Transparent		
3	1	Valve Membrane V2000 IV, 60A Durometer	Silicone, Transparent		
4	рс	Solvent	Tetrahydrofurane		

## UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE: VALUE PLASTICS, INC FRACTIONS DECIMAL ANGULAR † .XX ±.01 ± 1° .XXX ±.005 .XXXX ±.001 PART NUMBER: SCV67230/ML DESCRIPTION: Check Valve, Barb for 1/8" (3.2 mm) ID Tubing to Male Locking Luer DO NOT SCALE DRAWING DRAWN BY: DWG. NO. DATE: SIZE: SCALE: kyles 6/7/2013 B NTS 1027179 MANAGER: RaviN 1 OF 1 © 2011 VALUE PLASTICS, INC

## FOR REFERENCE ONLY