

NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE DIMENSIONS MARKED - ARE CONTROL DIMENSIONS


> ALL INTERIOR SURFACES HAVE A 15 Ra MICROINCH SURFACE FINISH OR FINER ALL EXTERIOR SURFACES HAVE A 20 Ra MICROINCH SURFACE FINISH OR FINER


REV.	DESCRIPTION
А	INITIAL RELEASE

	REVISIONS								
Ī	REV.	DESCRIPTION	DATE	APPROVED					
Ī	А	INITIAL RELEASE	01/16/12	S.N.R					









SECTION A-A

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE:		VALUE PLASTICS, INC						
FRACTIONS DECIMAL ANGULAR  ± .XX ±.01 ± 1° .XXX ±.005 .XXXX ±.001  DO NOT SCALE DRAWING		PART NUMBER: SFBP680						
		DESCRIPTION: SANITARY FITTING, BAG PORT 600 SERIES BARB, 1/2" (12.7mm) I.D. TUBING						
DRAWN BY:	DATE:	SIZE:	SCALE:	VER:	DWG. NO.			REV
kyles	1/17/2012	В	NTS		1022166			А
NEW PRODUCT DEVELOPMENT MGR:								
R. Narayanan	©	2010	VALUI	E PLASTICS, 1	INC SHEET:	1 OF 1		