




- NOTES: 1. FOR ADDITIONAL DETAILS, PLEASE CONTACT VP PRODUCT DEVELOPMENT  
 2. DIMENSIONS MARKED  $\text{---}$  ARE CONTROL DIMENSIONS  
 3. FERRULE END CAPS ARE DESIGNED TO WORK WITH COMMON COMPONENTS USED IN BIOPROCESSING INDUSTRY

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	INITIAL RELEASE	06-06-19	FJL



SECTION A-A



UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE:		<b>VALUE PLASTICS, INC</b>			
FRACTIONS    DECIMAL    ANGULAR ±            .XX ±.01            ± 1° .XXX ±.005 .XXXX ±.001		PART NUMBER: SFM4C			
DO NOT SCALE DRAWING		DESCRIPTION: Sanitary Fitting, Cap With 4" Flange			
DRAWN BY: Frank.Lombardi	DATE: 6/6/2019	SIZE: <b>B</b>	SCALE: NTS	VER: 	DWG. NO. 1082754
NEW PRODUCT DEVELOPMENT MGR: Doug Kimmel		© 2007 VALUE PLASTICS, INC			REV A
SHEET:				1 OF 1	

**FOR REFERENCE ONLY**