

A1	B1	C1	D1	E1	F1
<p>FCFML-2A-001 Stopcock, 1-way, Female Luer to Male Luer Body w/Positive stop, Clear PC Body, 2 Arrow Blue HDPE Diverter</p>	<p>FCHPFMR-002 High Pressure Stopcock, 500 psi, 1-way, Female Luer to Male Luer w/Luer Lock Ring, Clear PC Body, 2 Arrow Gold POM Diverter, Clear PC Lock Ring</p>	<p>FCFFM-014 Stopcock, 3-way, Female Luer to Female Luer to Male Luer, Rad. Stable PC Body, 1 Arrow Blue HDPE Diverter</p>	<p>FC4LB3F-001 Stopcock, 4-way, Female Luer to Female Luer to Elongated Female Luer, Large Bore, Rad. Stable PC Body, 1 Arrow White HDPE Diverter</p>	<p>FCHPFFMR-002 High Pressure Stopcock, 500 psi, 3-way, Female Luer to Female Luer to Male Luer w/Luer Lock Ring, Clear PC Body, 4 Arrow Gray Acetal Diverter, Clear PC Lock Ring</p>	<p>FCA026-001 Stopcock Assembly, 3-way, Female Luer to Female Luer to Bondable Socket w/7.2in Tubing Extension, Annealed Clear PC Body, 1 Arrow Natural HDPE Diverter, PVC Tube</p>
A2	B2	C2	D2	E2	F2
<p>FCFM-006 Stopcock, 1-way, Female Luer to Male Luer, Rad. Stable PC Body, 2 Arrow Natural HDPE Diverter</p>	<p>FCA019-001 High Flow Stopcock, 1-way, Female Luer to Male Luer w/Luer Lock Ring, Clear PC Body, 2 Arrow Natural HDPE Diverter, Rad. Stable PC Lock Ring</p>	<p>FCA018-001 Stopcock, 1-way, Female Luer to Threaded Post, Amber PEI Body, 2 Arrow Natural HDPE Diverter</p>	<p>FCFFM-1A-003 Stopcock, 3-way, Female Luer to Female Luer to Male Luer, Natural Nylon Body, 1 Arrow White HDPE Diverter</p>	<p>FCFFMR-004 Stopcock, 3-way, Female Luer to Female Luer to Male Luer w/Luer Lock Ring, Clear PC Body, 1 Arrow Natural HDPE Diverter, Clear PC Lock Ring</p>	<p>FCA047-001 Stopcock Assembly, 3-way, Female Luer to Female Luer to Bondable Socket, w/3in Tubing Extension and Bondable Male Luer, Annealed Rad. Stable PC Body, 1 Arrow Blue HDPE, PVC Tube, Annealed Clear PC Bondable Female Luer</p>
A3	B3	C3	D3	E3	F3
<p>FCFM-2A-004 Stopcock, 1-way, Female Luer to Male Luer, Natural Nylon Body, 2 Arrow Natural HDPE Diverter</p>	<p>FCFMR-001 Stopcock, 1-way, Female Luer to Male Luer w/Luer Lock Ring, Clear PC Body, 2 Arrow Natural HDPE Diverter, Clear PC Lock Ring</p>	<p>FCHPFFMR-005 High Pressure Stopcock, 500 psi, 1-way, Female Luer to Male Luer w/Luer Lock Ring, Clear PC Body, 2 Arrow Gray POM Diverter, Clear PC Lock Ring</p>	<p>FCFFM-002 Stopcock, 3-way, Female Luer to Female Luer to Bondable Socket, Clear PC Body, 1 Arrow Red HDPE Diverter</p>	<p>FCLBFFM-003 High Flow Stopcock, 3-way, Female Luer to Female Luer to Bondable Socket, Clear PC Body, 1 Arrow Natural HDPE Diverter</p>	<p>FCHPFFMR-003 High Pressure Stopcock, 808 psi, 3-way, Female Luer to Female Luer to Male Luer w/Luer Lock Ring, Clear PC Body, 4 Arrow Gold POM Diverter, Clear PC Lock Ring</p>