


| REVISIONS | | | |
|-----------|-------------|------|----------|
| REV | DESCRIPTION | DATE | APPROVED |
| | | | |



INSTRUMENT QUALITY PART.
NO PARTING LINE ALLOWED ON
THIS CRITICAL SURFACE.

| | | | |
|---|------------------------------------|---|-------------------------|
| UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE: | | VALUE PLASTICS, INC. | |
| FRACTIONS ± | DECIMALS .XX ±.01 .XXX ±.005 | ANGULAR ±1° | PART NUMBER: MLSL007 |
| DO NOT SCALE DRAWING | | DESCRIPTION: Male Luer to 500 Series Barb, 3/32" (2.25 mm) ID Tubing (May be used with separate stationary lock ring: FSLLR) | |
| DRAWN BY: SUE CORAM | SCALE: NTS | DRAWING NUMBER: | REV. |
| DATE: 08/21/97 | © 1997 VALUE PLASTICS, INC. | | SHEET: |

FOR REFERENCE ONLY