

NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE DIMENSIONS MARKED ARE CONTROL DIMENSIONS

	REVISIONS				
RE	v.	DESCRIPTION	DATE	APPROVED	
A	A	INITIAL RELEASE	06/07/13	S.N.R	





SCV67220/230 Bill Of Materials					
POS.	QTY	DESCRIPTION	MATERIAL, COLOR		
1	1	Check Valve V2000 Inlet with Barb 3/32" (2.4 mm) I.D. Tubing	MABS White		
2	1	Check Valve V2000 Outlet with Barb 1/8" (3.2 mm) I.D. Tubing	MABS, Transparent		
3	1	Valve Membrane V2000 IV, 60A Durometer	Silicone, Transparent		
4	рс	Solvent	Tetrahydrofurane		

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE: VALUE PLASTICS, INC FRACTIONS DECIMAL ANGULAR † .XX ±.01 ± 1° .XXX ±.005 .XXXX ±.001 PART NUMBER: SCV67220/230 DESCRIPTION: Check Valve, Barb 3/32" (2.4 mm) ID Tubing to Barb 1/8" (3.2 mm) ID Tubing DO NOT SCALE DRAWING DRAWN BY: SIZE: DWG. NO. DATE: SCALE: VER: kyles 6/7/2013 B NTS 1027176 MANAGER: SHEET: 1 OF 1 RaviN © 2011 VALUE PLASTICS, INC

FOR REFERENCE ONLY