

NOTES: FOR ADDITIONAL DETAILS, PLEASE CONTACT VP CUSTOMER SERVICE DIMENSIONS MARKED ARE CONTROL DIMENSIONS

REVISIONS								
REV.	DESCRIPTION	DATE	APPROVED					
А	INITIAL RELEASE	06/07/13	S.N.R					







SCV67FL/220 Bill Of Materials								
POS.	QTY	DESCRIPTION	MATERIAL, COLOR					
1	1	Check Valve V2000 Inlet with Female Luer Lock	SAN White					
2	1	Check Valve V2000 Outlet with Barb 3/32" (2.4 mm) I.D. Tubing	MABS, Transparent					
3	1	Valve Membrane V2000 IV, 60A Durometer	Silicone, Transparent					
4	рс	Solvent	Tetrahydrofurane					

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE: FRACTIONS DECIMAL ANGULAR ± .XX ±.01 ± 1° .XXX ±.005 .XXXX ±.001 DO NOT SCALE DRAWING		VALUE PLASTICS, INC								7		
		PART NUMBER: SCV67FL/220										
		Check Valve, Female Locking Luer to Barb for 3/32" (2.4 mm) ID Tubing										
DRAWN BY:	DATE:	SIZE:	SCALE:	VER:	DWG. NO							REV
kyles 6/7/2013		В	NTS	1027180				А				
MANAGER:												
RaviN		©	2011	VALUI	E PLAST	ICS,	INC	SHEET:	1	OF	1	

FOR REFERENCE ONLY