




NOTES: 1. FOR ADDITIONAL DETAILS, PLEASE CONTACT VP NEW PRODUCT DEVELOPMENT
2. DIMENSIONS MARKED  ARE CENTRAL DIMENSIONS
3. VENDOR CERTIFICATION OF COMPLIANCE REQUIRED
4. PART SHOWN MAY VARY SLIGHTLY IN APPEARANCE, MINOR DIFFERENCES NOT EFFECTING DIMENSIONAL CONFORMANCE ARE ACCEPTABLE

REVISIONS			
REV.	DESCRIPTION	DATE	APPROVED
A	INITIAL RELEASE	10/02/14	H.L
B	ECO 345	03-27-17	FJL
C	ECO 384	06-20-18	FJL



ASSEMBLY SKU	ITEM	DESCRIPTION	QTY
TBA2025-003	1	ADAPTER BODY W/OUT SIDEPORT - DEHP FREE, CLEAR PVC	1
	2	GASKET .025 - NATURAL SILICONE	1
	3	MALE CAP - BLUE PA ABS	1
TBA2025-004	1	ADAPTER BODY W/OUT SIDEPORT - DEHP FREE, CLEAR PVC	1
	2	GASKET .025 - SILICONE	1
	3	MALE CAP - BLUE PA ABS	1

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ARE: FRACTIONS DECIMAL ANGULAR ± .XX ±.01 ± 1° .XXX ±.005 .XXXX ±.001 DO NOT SCALE DRAWING		VALUE PLASTICS, INC						
		PART NUMBER: TBA2025						
		DESCRIPTION: Tuohy Borst Adapter without side port, Body, Cap, 0.025" ID Gasket, DEHP-Free						
DRAWN BY: marie.morleyjahn		DATE: 6/20/2018		SIZE: B	SCALE: NTS	VER: 	DWG. NO. 1032903	REV C
MANAGER: Frank.Lombardi		© 2017 VALUE PLASTICS, INC						SHEET: 1 OF 1

FOR REFERENCE ONLY