

Compliant to the ISO 80369-3 Standard

ENT Series connectors were designed, tested, qualified, and documented to be compliant with the recently published ISO 80369 standard. Small, lightweight, and easy to clean, these connectors offer comfort for everyone from neonatal infants to the elderly. The medical-grade materials are clear with no colorants added, and are lipid-resistant. Available in a wide range of sizes, ENT Series connectors can also be customized for additional non-standard sizes and colors.



Features & Benefits

- ISO 80369-3 or ENFit™ compliant
- Small, lightweight, easy to clean design
- One-piece tethered caps and plugs
- Lipid-resistant materials
- Clear materials with no colorants
- Extension assemblies and custom sizes available

Specifications

- **Materials:**
Medical-grade natural ABS
Also available in orange and purple
- **Sizes:**
10 FR female connectors
5 FR, 6.5 FR, 8 FR, & 10 FR male sizes for bond-in tubing
- **Barb Sizes:**
1/8" (3.2 mm) ID barb for silicone tubing

ENT Series

Enteral Feeding Connectors



ENTME050 (-8003)

Enteral (ENFit™) Male Connector with 5 FR Tubing Bond-in Pocket in Natural ABS



ENTME098 (-8003)

Enteral (ENFit™) Male Connector with 10 FR Tubing Bond-in Pocket in Natural ABS



ENTME064 (-8003)

Enteral (ENFit™) Male Connector with 6.5 FR Tubing Bond-in Pocket in Natural ABS



ENTM230 (-6005)

Enteral (ENFit™) Male Connector Barb for 1/8" (3.2 mm) ID Tubing in Animal-Free Polypropylene



ENTME078 (-8003)

Enteral (ENFit™) Male Connector with 8 FR Tubing Bond-in Pocket in Natural ABS



ENTFE098 (-8003)

Enteral (ENFit™) Female Connector with 10 FR Tubing Bond-in Pocket in Natural ABS

ISO 80369-3 Accessories

Tethered caps and plugs are not covered under the ISO 80369-3 standard, but are compatible with our enteral male and female connectors.



ENTFTC (-VP1)

Tethered Cap Fits Enteral (ENFit™) Male Connector in Animal-Free Polyethylene



ENTMTP (-VP1)

Tethered Plug Fits Enteral (ENFit™) Female Connector in Animal-Free Polyethylene

MAR-ENT-DS-01

Nordson MEDICAL
805 West 71st Street
Loveland, CO 80538, U.S.A.

+1.888.404.5837
+1.970.267.5200
nordsonmedical.com

Nordson
MEDICAL